



Reagent powder pillows, molybdenum 1, 50 mL, pk/100

Product #: 2352769
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed package for the determination of Molybdenum, 0.02-3.00mg/l Mo.

For determination of low range molybdenum by the ternary complex method (Hach method 8169). Range: 0.02-3.00 mg/L Mo Sample Size: 50 mL

Specifications

Method:	Ternary Complex
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)