



MolyVer 2 Molybdenum Reagent Powder Pillows, pk/100

Product #:	2604399
ZAR Price:	Contact Hach
Ships within 3 days	

Powder reagent in sealed package for the determination of Molybdenum, 0.3-40.0mg/L Mo. For 10mL sample.

Specifications

Method:	Mercaptoacetic Acid
Number of tests:	100
Parameter:	Molybdenum
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)