



Reagent set, LR molybdenum, 5 mL

Product #:	2439500
ZAR Price:	Contact Hach
Available	

For determination of Low Range Molybdate Molybdenum by the ternary complex method, using visual test kits. Replacement reagent set for MO-LR Test Kit (2359300). Range: 0.25 - 3 mg/L Mo. Sample Size: 5 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Specifications

Method:	Ternary Complex
Range:	0.25 - 3 mg/L Mo
Sample size:	5 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)