



Reagent set, molybdenum, LR, 50 mL

Product #: 2440200
ZAR Price: Contact Hach
Ships within 3 days

For determination of Low Range Molybdate Molybdenum by the Ternary Complex method (Hach method 8169). Range: 0.02 - 3.00 mg/L as Mo. Sample Size: 50 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Specifications

Instrument:	DR/3000, DR/2000
Method:	Ternary Complex
Number of tests:	100
Range:	0.02 - 3.00 mg/L Mo
Sample size:	50 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)