



Reagent solution, molybdenum 2, 50 mL MDB

Product #: 2352512
ZAR Price: Contact Hach
Ships within 3 days

For determination of low range molybdenum by the ternary complex method (Hach method 8169). Range: 0.02-3.00 mg/L Mo 50 mL Marked Dropping Bottle

Specifications

Method:	Ternary Complex
Number of tests:	100
Platform:	Solution
Quantity:	50

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)