



Molybdenum 2 Reagent Solution (100mL MDB)

| | |
|---------------------|--------------|
| Product #: | 2352532 |
| ZAR Price: | Contact Hach |
| Ships within 3 days | |

For the determination of Molybdenum, 0.02-3.00mg/L Mo. For 50mL sample. Measured dropping bottle.

Specifications

| | |
|------------------|-----------------|
| Method: | Ternary Complex |
| Number of tests: | Varies |
| Platform: | Powder Pillows |
| Quantity: | 100 mL |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)