



Molybdenum Reagent Set, HR, Powder Pillows, 0.2 - 40.0 mg/L Mo

Product #: 2604100
ZAR Price: Contact Hach
Ships within 3 days

Powder reagents in sealed packages for the determination of Molybdenum by the Mercaptoacetic Acid Method. For 10mL sample. Pack of 100.

Specifications

Instrument:	All except DR/3000, DR/2000, PC II
Method:	Mercaptoacetic Acid
Number of tests:	100
Parameter:	Molybdenum
Platform:	Powder Pillows
Range:	0.3 - 40.0 mg/L Mo
Sample size:	10 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)