



Nickel Reagent Set, Heptoxime Method Nickel Reagent Set, 0.02-1.80mg/L Ni

Product #:

2243500

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For the determination of Nickel, Method: Chloroform extraction using Heptoxime Indicator.

Specifications

EPA compliant:	Yes
Instrument:	All except DR 900, DR/800 Series, PC II
Method:	Heptoxime
Number of tests:	50
Parameter:	Nickel
Range:	0.02 - 1.80 mg/L Ni

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)