



Be Right™



Potassium 1 Reagent Powder Pillows, pk/25

Product #: 1432198
ZAR Price: Contact Hach
Available

Powder reagent in sealed packages for the analysis of Potassium, 0.1-7.0mg/L K. For 25mL sample.

Specifications

Method: 8049
Number of tests: 25
Parameter: Potassium
Platform: Powder Pillows
Quantity: 25 /pk

Required Accessories

- Potassium 3 Reagent Powder Pillows, pk/100 (Item 1432399)
- Potassium 2 Reagent Solution Pillows, pk/25 (Item 1432298)
- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)