



## Potassium 3 Reagent Powder Pillows, pk/100

**Product #:**

**1432399**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagent in sealed package for the analysis of Potassium, 0.1-7.0mg/l K.

For potassium determination in water, soil and plant tissues by the tetraphenylborate method. Hach Method 8049 Sample size: 25 mL. Pack of 100 powder pillows.

---

### Specifications

Method:	Tetraphenylborate
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)