



Potassium Reagent Set, 0.1-7.0mg/L K

Product #:

2459100

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Powder reagents and reagent solution in sealed packages for determination of Potassium by the Tetraphenylborate method. For 25mL samples; pk/100.

Specifications

Description:	Potassium Reagent Set
Instrument:	All except DR800 Series, PC II
Method:	Tetraphenylborate
Number of tests:	100
Parameter:	Potassium
Platform:	Powder Pillows
Range:	0.1 - 7.0 mg/L K
Special Remarks:	Set contains: Potassium 1 Reagent Powder Pillows (1432198), Potassium 2 Reagent Solution Pillows (1432298), Potassium 3 Reagent Powder Pillows (1432399)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)