



Reagent powder pillows, potassium 3, pk/1000

Product #: 1432328

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:

- 1432399

If necessary, we will change your selected shipping method to accomodate these items.

For potassium determination in water, soil and plant tissues by the tetraphenylborate method. Hach Method 8049 Sample size: 25 mL. Pack of 1000 powder pillows.

Specifications

Method:	Tetraphenylborate
Platform:	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)