



Potassium 2 Reagent Solution Pillows, pk/25

Product #:

1432298

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Solution pillows in sealed packages for the analysis of Potassium, 0.1-7.0 mg/L K. For 25mL sample.

Specifications

Method:	Tetraphenylborate
Number of tests:	25
Parameter:	Parameter
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)