



Silver Reagent Set, Powder Pillows, 0.02-0.700 mg/L Ag

Product #: 2296600
ZAR Price: Contact Hach
Ships within 3 days

Powder reagents in sealed packages for colorimetric Silver Method using Cadion 2B indicator. For 50mL samples.

Specifications

Instrument:	All except DR 900, DR/800 Series, PC II
Method:	8120
Number of tests:	50
Parameter:	Silver
Platform:	Powder Pillows
Range:	0.02 - 0.70 mg/L Ag

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)