



Reagent powder pillows, silver 1, pk/50

Product #:	2293566
ZAR Price:	Contact Hach
Ships within 3 days	

Reagent for determination of Silver by the Cadion 2B Colorimetric Method. Hach Method 8120. Range: 0.02 - 0.70 mg/L Ag. Sample Size: 50 mL. Pack of 50 Pillows. Also available in reagent set 2296600.

Specifications

Method:	Colorimetric
Number of tests:	50
Platform:	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)