



Reagent solution pillows, silver 2, pk/50

Product #: 2293666
ZAR Price: Contact Hach
Ships within 3 days

Reagent for determination of Silver by the Cadion 2B Colorimetric Method. Hach Method 8120. Range: 0.02 - 0.70 mg/L Ag. Sample Size: 50 mL. Pack of 50 Pillows. Also available in reagent set 2296600.

Specifications

Method: Colorimetric
Number of tests: 50

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)