



## Zinc Reagent Set, 0.01-2.00mg/L Zn

Product #:

2429300

ZAR Price:

Contact Hach

Ships within 3 days



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Zinc, Zincon Method. For 20 mL sample. pk/100. (Digestion is required for determining total Zinc). Sample Size: 20 mL. Range: 0.01 - 3.00 mg/L Zn.

---

### Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Method:	8009
Number of tests:	100
Parameter:	Zinc
Range:	0.01 - 3.00 mg/L Zn

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)