



Ammonia Salicylate Reagent Powder Pillows

Product #:	2653299
ZAR Price:	Contact Hach
Ships within 3 days	

Powder reagent in sealed packages for the determination of Nitrogen, Ammonia (Salicylate Method), 0.01-0.50mg/L $\text{NH}_3\text{-N}$. For 10 mL sample. 100/pk.

Specifications

Method:	8155
Number of tests:	100
Parameter:	Ammonia
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	0.01 - 0.50 mg/L $\text{NH}_3\text{-N}$