



Aqua Ammonia Reagent Set, Digital Titrator

Product #: 2930500

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Set for Model NI-HRDT Digital Titrator Kit, Cat. No. 2930400. Range: 5 - 35% (50 - 350 g/L) NH_3 . Includes Sulfuric Acid Digital Titrator Cartridge and Wide Range pH Indicator. Approximately 100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	10222
Number of tests:	100
Platform:	Digital Titrator
Range:	50 - 350 g/L NH_3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)