



## Nitrogen, Ammonia TNT, 0.02 - 2.50 mg/L NH<sub>3</sub>-N

Product #:

2604545

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests

---

### Specifications

|                     |  |
|---------------------|--|
| Method:             | Salicylate   |
| Number of tests:    | 25 - 50  |
| Parameter:          | Ammonia  |
| Platform:           | Test 'N Tube   |
| Range:              | 0.02 - 2.50 mg/L NH <sub>3</sub> -N  |
| Special Remarks:    | Set contains: Cuvette test tubes, AmVer Diluent Reagent (2602200), Ammoniasalicylate Powder Pillows (2395266), Ammoniacyanurate Powder Pillows (2395466) |
| Storage conditions: | 10 - 25 °C   |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)