



Nitrogen, Ammonia TNT, 0.4 - 50 mg/L $\text{NH}_3\text{-N}$

Product #:

2606945

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube 25-50 tests.

Specifications

Method:	10031
Number of tests:	25 - 50
Parameter:	Ammonia
Platform:	Test 'N Tube
Range:	0.4 - 50.0 mg/L $\text{NH}_3\text{-N}$
Special Remarks:	Set contains: Cuvette test tubes, AmVer Diluent reagent, Ammoniasalicylate Powder Pillows (2395266), Ammoniacyanurate Powder Pillows (2395466)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)