



Nitrate Reagent Set, Powder Pillows, 0.01-0.50 mg/L NO₃-N

Product #: 2429800
ZAR Price: Contact Hach
Ships within 3 days

Powder reagents in sealed packages for the determination of Nitrogen, Nitrate by the Cadmium reduction method. 10mL sample, 100 pcs.

Specifications

Instrument:	All except DR/3000, DR/2010, DR/2000, PC II
Method:	Cadmium Reduction
Number of tests:	100
Parameter:	Nitrate, Nitrogen
Platform:	Powder Pillows
Range:	0.01 - 0.50 mg/L NO ₃ -N
Special Features:	Nitrate Reagent set, Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)