



Nitrate Reagent Powder Pillows, 0.3-30mg/L NO₃-N

Product #:

2106169

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Powder reagent in sealed packages for the determination of Nitrate by Cadmium reduction method. For 10 mL sample, 100 pieces.

Specifications

Description:	NitraVer 5 Nitrate Reagent Powder Pillows, 10 mLHR, MR
Instrument:	All except DR/3000, DR/2010, DR/2000
Method:	8039 HR
Method Number 2:	8171 MR
Number of tests:	100
Parameter:	Nitrate
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	0.3 - 30.0 mg/L NO ₃ -N
Range 2:	0.1 - 10.0 mg/L NO ₃ -N
Reference method:	SM 4500-NO ₃ -B
Storage conditions:	10 °C to 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)