



Nitrate TNT, 0.2 - 30 mg/L NO₃-N

Product #:

2605345

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Method:	10020
Number of tests:	50
Parameter:	Nitrate
Platform:	Test 'N Tube
Quantity:	50 /pk
Range:	0.2 - 30.0 mg/L NO ₃ -N
Special Remarks:	Set contains: Cuvette test tubes, NitraVer X Reagent A (2605400), NitraVer X Reagent B Powder Pillows (2605546)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)