



Nitrite Reagent Powder Pillows, 2-250 mg/L NO₂

Product #: 2107569
ZAR Price: Contact Hach
Ships within 3 days

Powder reagents in sealed packages for the determination of Nitrite by the Ferrous Sulfate Method. For 10mL sample. Pack of 100.

Specifications

Instrument:	All except DR/3000, DR/2000, DR/850, DR/820, PC II
Method:	8153
Number of tests:	100
Parameter:	Nitrite
Platform:	Powder Pillows
Quantity:	100 /pk
Range:	2 - 250 mg/L NO ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)