



NitriVer 3 Nitrite Reagent Powder Pillows, 5 mL, pk/100

Product #:	1407899
ZAR Price:	Contact Hach
Ships within 3 days	

For nitrite determination by the NitriVer 3 diazotization method. Sample size 5 mL. For use with Hach visual test kits. Range: up to 0.5 mg/L NO₂-N.

Specifications

Method:	Diazotization
Platform:	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)