



Nitrogen, total, TNT, 0.5-25 mg/L N

Product #:

2672245

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Digestion Required:	Yes
Method:	Persulfate Digestion
Number of tests:	25 - 50
Parameter:	Nitrogen, total
Platform:	Test 'N Tube
Range:	0.5 - 25.0 mg/L N
Special Remarks:	Set consists of 2 cartons: 26717-45 and 26721-45 26717-45 = TN Hydroxide reagent (2714045), Persulfate Powder Pillows (2671846) 26721-45 = TN Acid reagent C (2672145), Bisulfite reagent A Powder Pillows (2671946), Indicator reagent B Powder Pillows (2672046)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)