



Nitrogen, total, TNT, 10-150 mg/L N

Product #:

2714100

 **Hazardous**

ZAR Price:

Contact Hach

Ships within 3 days

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests

Specifications

Digestion Required:	Yes
Method:	Persulfate Digestion
Number of tests:	25 - 50
Parameter:	Nitrogen, total
Platform:	Test 'N Tube
Range:	2 - 150 mg/L N
Special Remarks:	Set consists of 2 cartons: 27140-45 and 26721-45 27140-45 = TN Hydroxyde reagent solution (2714045), Persulfate Powder Pillows (2671846) 26721-45 = TN Acid reagent C solution (2672145), Bisulfite reagent A Powder Pillows (2671946), Indicator reagent B Powder Pillows (2672046)
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)