



Sodium Hydroxide Digital Titrator Cartridge, 0.3636 N

Product #:

1437801

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	8205
Number of tests:	40 - 100
Parameter:	Carbon Dioxide
Platform:	Digital Titrator
Range:	10 - 100 mg/L CO ₂
Titrant:	0.3636 N NaOH

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)