



Anionic Surfactants TNTplus Vial Test (0.1-4.0 mg/L), 25 Tests

Product #:
ZAR Price:
No availability shown

TNT874
Contact Hach



Hazardous
Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Expert Anionic Surfactants test made simple.

Anionic Surfactants TNTplus chemistry and your Hach spectrophotometer are engineered to simplify water analysis for accurate results, everytime.

For determination of Anionic Surfactants by the Methylene Blue method.

Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

Specifications

Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Method:	10278
Number of tests:	25
Parameter:	Surfactants, Anionic
Platform:	TNTplus
Range:	0.1 - 4.0 mg/L
Standard method:	ISO 7875-1-2-1984, DIN 38409-H 23-2
What's included?:	25 test vials (3 mL), 1 reagent bottle (10 mL)

What's included?

25 test vials (3 mL), 1 reagent bottle (10 mL)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)