



Phenol 2 Reagent Powder Pillows

Product #: 183699
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed packages for the determination of phenols by the 4-Aminoantipyrine method. pk/100

Specifications

Description:	Phenol 2 Reagent Powder Pillows
Method:	8047
Number of tests:	50
Parameter:	Phenols
Platform:	Powder Pillows
Quantity:	100
Storage conditions:	10 °C to 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)