



## Phenols Reagent Set, 0.002-0.200 mg/L

**Product #:**

**2243900**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Phenols, 4-Aminoantipyrine Method.

---

### Specifications

Description:	Phenols Reagent Set
EPA compliant:	Yes
Instrument:	All except DR 900,  DR/800 Series, PC II
Method:	4-Aminoantipyrine
Number of tests:	100
Parameter:	Phenols
Range:	0.002 - 0.200 mg/L
Shipping Method:	No Air Shipments

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)