



Phenols Reagent Set (Air Shipment)

Product #:

2243901

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Phenols by the 4-Aminoantipyrine method. Hach Method 8047. USEPA accepted for reporting wastewater analyses. Range: 0.002 - 0.200 mg/L Reagent Set includes: Hardness 1 Buffer Solution, Chloroform, Phenol Reagent Powder Pillows, and Phenol 2 Reagent Powder Pillows for 100 tests. Reagents packaged for air freight shipments.

Specifications

| | |
|------------------|--|
| Description: | Phenols Reagent Set |
| EPA compliant: | Yes |
| Instrument: | All except DR 900, DR/800 Series, PC II |
| Method: | 4-Aminoantipyrine |
| Number of tests: | 100 |
| Range: | 0.002 - 0.200 mg/L |
| Shipping Method: | Air Shipments |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)