



QAC reagent 1, 0.2- 5.0 mg/L CTAB, 50 powder pillows

Product #: 2401066
ZAR Price: Contact Hach
Ships within 3 days

Powder reagent in sealed package for the determination of Quaternary Ammonium Compounds (QAC).

Specifications

Method:	8336
Number of tests:	50
Parameter:	Quaternary Ammonium Compounds
Platform:	Powder Pillows
Quantity:	50/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)