



Tannin & Lignin Reagent Set, 0.1-9.0 mg/L

Product #:

2244600

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Tannin and Lignin by the Tyrosine Method. For 25 mL sample.

Specifications

Instrument:	All except DR 900, PC II
Method:	8193
Number of tests:	100
Parameter:	Tannin & Lignin
Platform:	Solution
Range:	0.1 - 9.0 mg/L as Tannic Acid

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)