



## Reagent, TanniVer 3, tannin-lignin, 100 mL MDB

**Product #:**

**256032**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of tannin & lignin by the Tyrosine method

---

### Specifications

Method:	Tyrosine
Number of tests:	100
Platform:	Solution

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)