



Standard solution, sulfuric acid, 19.2 N, 500 mL

Product #:

203849

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Multi-purpose sulfuric acid standard solution. Required for determination of volatile acids by the esterification method. 500 mL bottle

Specifications

Concentration:	19.2 N
Method:	Esterification
Platform:	Solution
Volume:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)