



## Nitrification inhibitor for BOD (35 g)

**Product #:** 253335  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

TCMP - 2-chloro-6 (trichloromethyl) pyridine (N-Serve), coated on an inert substrate. Use 0.16 gr in a 300-mL sample inhibits nitrogenous oxygen demand in the BOD test (Use Dispenser Cap, Part No. 45901), 35 g.

---

### Specifications

Application:	Nitrification Inhibitor  (CBOD)
Description:	Nitrification Inhibitor
Formula:	Trichloromethyl Pyridine (TCMP)
Method:	8043
Number of tests:	200
Package Type:	Bottle
Parameter:	Biochemical Oxygen Demand
Platform:	Jar
Quantity:	35 g

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)