



BOD Nutrient Buffer Pillows, 4 mL concentrate to prepare 4 L, 50/pk

Product #:

2436466

ZAR Price:

Contact Hach

Ships within 3 days

BOD Reagents

For preparation of BOD5 dilution water. Each pillow combines the APHA-specified buffer and nutrients to prepare 4 liters of dilution water. Pack of 50 pillows.

Specifications

Application:	Dilution water preparation
Description:	Nutrients & Buffer
Formula:	Nutrients & Buffer
Package Type:	Powder Pillows
Quantity:	4 mL
Units:	50/pk

Required Accessories

- BOD Standard Solution, 300 mg/L, pk/16 - 10-mL Voluette Ampules (Item 1486510)
- Nitrification inhibitor for BOD (35 g) (Item 253335)
- BOD Seed Inoculum, Polyseed, pk/50 capsules (Item 2918700)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)