



## Buffer solution pillows, respirometric BOD nutrient, pk/50

**Product #:** 2962266  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For use with BODTrak II Apparatus. Meets Standard Methods 5210 D requirements. Use one pillow per sample bottle. Included with BODTrak II Apparatus (2952400) and BODTrak II Reagent Set with Bottles (2952440).

---

### Specifications

Application:	BODTrak II
Description:	Nutrients & Buffer
Formula:	Nutrients & Buffer
Package Type:	Powder Pillows
Quantity:	50 pcs
	1 Pillow/Bottle

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)