



Filling Solution, Internal, Polarographic DO, 50 mL

Product #: 2759126
ZAR Price: Contact Hach
Available

Filling solution for Polarographic Dissolved Oxygen probe, Potassium Chloride/Glycerine, 50 mL bottle

Specifications

Description:	5130, 5197000
Formula:	KCl & Glycerine
Platform:	Solution
Quantity:	50 mL
Solution Type:	Refill

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)