



Dissolved Oxygen Reagent Accuvac Ampules, 6-800 µg/L O₂

Product #: 2501025
ZAR Price: Contact Hach
Ships within 3 days

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Dissolved Oxygen, Indigo Carmine Method.

Specifications

Instrument:	All except DR 820, PC II
Method:	8316
Number of tests:	25
Parameter:	Dissolved Oxygen
Platform:	Accuvac
Quantity:	25 /pk
Range:	6 - 800 µg/L O ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)