



Dissolved Oxygen Reagent Accuvac Ampules, 0.3-15.0 mg/L O₂

Product #: 2515025
ZAR Price: Contact Hach
Ships within 3 days

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Dissolved Oxygen, HRDO Method.

Specifications

Instrument:	0.3-15.0 mg/L: All 1.0-40.0 mg/L: All except DR 900, DR/800 series
Method:	HRDO
Method Number 2:	8166, 8333
Number of tests:	25
Parameter:	Dissolved Oxygen
Platform:	Accuvac
Quantity:	25 /pk
Range:	0.3 - 15.0 mg/L O ₂
Range 2:	1.0 - 40.0 mg/L O ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)