



Water, deionised and demineralised, 500mL

Product #: 27249
ZAR Price: Contact Hach
Ships within 3 days

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specifications

Description:	Water, Deionized
Platform:	Solution
Volume:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)