



Formazin Turbidity Standard, 4000 NTU, 100 mL

Product #: 246142
ZAR Price: Contact Hach
Ships within 3 days

Turbidity standard can be diluted to be used with any of Hach's turbidimeters TL23, TU52, 2100Q, 2100N, 2100AN series

Specifications

Concentration:	4000 NTU
Parameter:	Turbidity
Platform:	Solution
Quantity:	100 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)