



Stablcal Turbidity Standards Calibration Kit for 2100P

Product #:	2659405
ZAR Price:	Contact Hach
Available	

Four ready to use, stabilised Formazin Primary standard solutions. Kit includes sealed vials of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU.

Specifications

Concentration:	< 0.1, 20, 100, 800 NTU
Instrument:	2100P
Parameter:	Turbidity
Quantity:	Set of 4 sealed vials

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)