



Calibration set, Stablcal turbidity standards, 2100P portable turbidimeter, 100 mL bottles

Product #:	2659410
ZAR Price:	Contact Hach
Available	

Four ready to use, stabilised Formazine Primary standard solutions in bottles. Kit includes 100 mL of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

Specifications

Concentration:	< 0.1, 20, 100, 800 NTU
Quantity:	Set of 4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)