



Calibration kit, Stablcal turbidity standards, 2100AN / AN IS turbidimeter, 100 mL bottles

Product #: 2659510
ZAR Price: Contact Hach
Ships within 2 weeks

Ready for use, stabilised Formazine Primary Standards in bottles, Set: <0.1, 20, 200, 1000, 4000, 7500 NTU

Specifications

Concentration:	< 0.1, 20, 200, 1000, 4000, 7500 NTU
Description:	Stablcal Calibration Kit for 2100 AN / AN IS (100 mL Bottles)
Instrument:	2100 AN / AN IS
Parameter:	Turbidity
Platform:	Solution
Quantity:	Set of 6
Storage conditions:	5 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)