



Stablcal Stabilised Formazine Turbidity Standard #0.1 NTU, 500 mL

Product #:	2659749
ZAR Price:	Contact Hach
Ships within 3 days	

Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, <0.1 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Specifications

Concentration:	< 0.1 NTU
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)