



## Stablcal Stabilised Formazine Turbidity Standard, 20 NTU, 500 mL

<b>Product #:</b>	<b>2660149</b>
ZAR Price:	Contact Hach
Available	

Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, 20 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

### Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### Based on Formazin Primary Standards

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Specifications

Concentration:	20 NTU
Platform:	Solution

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)